

Trigger Mode Matlab

%% Set parameters

SetSingleChanAllParam_v2(s, Channel ID, t1, t2, t3, t4, N, I, 1) //

%% Set parameters

```
SetSingleChanAllParam_v2(s, 0, ...  
    pulse_width, ...    % pulseDurationUS  
    pulse_deadtime, ...    % deadTimeUS  
    interpulse_duration, ...    % interpulseDurationUS  
    interframe_duration, ...    % interframeDurationUS  
    N_pulse_repetition, ...    % numberOfPulsesPerFrame  
    current, ...    % IAmplitude in mA  
    0); %Trigger mode
```

With Continuous Trigger enabled (SetSingleChanSingleParam_v2(s, 0, 7, 0))

Before command for trigger: SetSingleChanState(s, 0, 1, 1, 1);

- **Command : SetSingleChanSingleParam_v2(s, 0, 9, 1):**
When executing this command, in which the stimulator is in a mode that operates in a context of continuously sent pulses (following the interframe time), the behavior is as follows: immediately after the execution of the command, a pulse will be retransmitted and **the interframe time interval will be reset**. Then, the **interframe time interval will be respected**.

With Continuous Trigger disabled (SetSingleChanSingleParam_v2(s, 0, 7, 1)).

- **Command: SetSingleChanSingleParam_v2(s, 0, 8, 1):**
In this scenario, the electrostimulator is in the mode where continuous pulses are not being generated. Thus, **each command sent 1 pulse**. Sending 2 commands results in two pulses being sent, **respecting the interframe time**.

Example: If 2 commands are sent within a period shorter than the interframe time, the second pulse will be generated only after the interframe time.

CAUTION: If the second command is not sent within the interframe time, the second pulse will not be generated.

Optional - Sending triggers without respecting time:

Command: SetSingleChanSingleParam_v2(s, 0, 8, 1)

Sends the trigger whenever a command is sent, without respecting the time between pulses. (Check if there is any usefulness)